



### Intended use

PULIH is intended to be used for monitoring rehabilitation progress of patients with chronic low-back pain under supervision of health-care providers, physicians, rehabilitation specialists, physiotherapist, chiropractors and normal user. PULIH will be used with PULIH mobile apps (iOS & Android) for control of the measurements and personal progress tracking. Clinicians will be able to monitor the patients progress remotely using the Dr Dashboard which can be accessed through internet browser. It is suitable to be used at home or healthcare facility. Before using the device, please review the Instruction for Use by scanning the QR code given.

**1**

electrode1  
electrode2  
electrode3

Slider housing  
Attach the electrode to the slider

**2**

L4  
L5

- Follow the recommended instruction by physician and place the slider at L4 and L5 position at the back, use the same position every time for consistent measurement result

**3**

- Put the device on the flat surface to start the connection
- Turn on the device power button, LED in blue color
- LED will blink for Bluetooth connection from Mobile Apps

**4**

- Launch the mobile apps, then the device will connect via Bluetooth.
- The blue LED light will stop blinking and be still light

**5**

- Click the "+" sign to start, then click "continue". Select "start" to activate for measurement.
- The device LED will turn into Green color and ready for measurement

**6**

body

- The device attach it to the slider in the correct magnetic polarity

**7**

- Follow the screen instruction and start the bowing

**8**

You will see the result once the measurement complete



**Package Specification:**  
Device, Slider, Charging Cable & Quick reference instruction



For Inquirees, please contact us at

✉ [info@nonivasi.com](mailto:info@nonivasi.com)

🌐 [www.nonivasi.com](http://www.nonivasi.com)

M1-08-01, 8 Trium Menara 1, Jalan Cempaka Sd 12/5 Bandar Sri Damansara, 52200 Kuala Lumpur, Malaysia